SCHILLER

MEDICAL S.A.S.

FRED easy TRAINER

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Contents

1	General points	20
2	Intended use	22
3	Characteristics	22
4	Display and controls	23
4.1	Device	23
4.2	Electrodes and cable	24
5	Scenario selection procedure	25
5.1	List of scenarios	25
5.2	Selection of a scenario.....	25
5.2.1	Selection of a scenario on the device.....	25
5.2.2	Selection of a scenario with the infrared remote control.....	26
5.3	Description of the functions associated with keys on the infrared remote control.....	27
5.4	To switch the device off	27
6	Safety instructions	28
7	Cleaning and maintenance	28
8	Servicing and disposal at the end of the service life	28
9	Specifications	29
9.1	Physical and environment requirements	30
10	Warranty	30
11	Accessories	31
12	General algorithm	32
13	Spoken and displayed texts	33

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1 General points

The FRED easy TRAINER device bears the CE mark in accordance with Directive 93/42/EEC relating to medical devices, and fulfils the fundamental requirements of Appendix I of that Directive.

- The device complies with standard EN 55011 class B.
- The safety of the user and trouble free device operation, shall only be guaranteed when original SCHILLER parts are correctly used. Further, only the accessories mentioned in these instructions for use or accessories that have been expressly authorised by SCHILLER may be used. If the device becomes damp as a result of temperature variations, it may only be put back into service once it is completely dry.
- The manufacturer shall not be held liable for the safety, reliability and performance of the device if:
 - Assembly, extensions, adjustments, modifications or repairs are not performed by the manufacturer or by personnel authorised by the manufacturer.
 - The electrical installation of the premises of use does not comply with locally applicable regulations.
 - The device is not used in accordance with the instructions for use.
- The manufacturer may supply, on request, circuit diagrams, lists of components, descriptions, calibration instructions or any other information, for use by the qualified technical personnel of the user for repairing such parts of the device as have been stated to be repairable by the manufacturer. Such supply shall not in any case constitute permission or approval to modify or repair the device.
- This document describes a device type that complies with the safety standards applicable at the time of printing. SCHILLER reserves all rights to the circuits, processes, names, software and devices mentioned in these instructions for use.
- The user guide is an integral part of the device and should always be kept near the device. Close observance of the information given in the user guide is a prerequisite for using the device as intended for and correct operation and ensures operator safety. Therefore, be sure to read the complete user guide.
- The user guide informs the device operator about the intended use, exact function, operation and required preventive maintenance. It is not a substitute for a product training.
- The quality assurance process used by SCHILLER is in compliance with international standards ISO 9001 and ISO 13485.
- No part of these instructions for use may be reproduced without the written permission of SCHILLER.
- This device is not designed for any use other than that specifically described in this manual. Such use may be hazardous

Manufacturer

SCHILLER Médical SAS
4, rue Louis Pasteur – BP 90050
F-67162 WISSEMBOURG Cedex
Tel: +33 (0)3.88.63.36.00
Fax: +33 (0)3.88.94.12.82
E-mail: info@schiller.fr
Internet site: www.schiller-medical.com

Revision history:		
Version	Date	Comment
01-01	24.04.2003	1st issue
Août 2004	August 2004	1st revision
September 2004	September 2004	addition of German version
March 2005	March 2005	addition of alkaline battery
April 2005	April 2005	modification of alkaline battery
May 2005	May 2005	Correction on side 20
May 2006	May 2006	New guidelines ERC / AHA 2005

2 Intended use

The FRED easy TRAINER is a simple cost-effective simulator of the FRED easy that is designed to be used exclusively for training. In no event may it be used as a conventional defibrillator. It should not be connected to a patient. It should only be used on a mannequin. It has been designed for the purpose of training users of SCHILLER's FRED easy automated defibrillator. The device provides several simulations or scenarios (from 0 to 8) that aim to help to familiarise with the FRED easy and shows that the users have the knowledge required to use it in case of emergency. The device also allows the trainer to work using an infrared remote control to select a scenario or to intervene in real time on the course of an intervention.

3 Characteristics

The FRED easy TRAINER allows the **simulated** delivery of electric shocks. There is no high voltage in the device, and training is risk free.

The FRED easy TRAINER may be used on all types of mannequins as it does not deliver an electric shock. It may only be used with the training electrodes sold by Schiller.

The FRED easy TRAINER may be used with an infrared remote control. The remote control allows the trainer to control the training actions.


The FRED easy TRAINER has 9 pre-recorded training scenarios which simulate realistic cardiac arrest configurations.

The basic equipment of the FRED easy TRAINER consists of :

- 1 rechargeable 6 Volts battery and a battery charger or 1 alkaline 6 Volts battery
- Training electrodes consisting of a single use set of 2 electrodes and a reusable cable. The connector of the cable is fitted with a magnet allowing the device to recognise the presence of the electrodes in the device without any electrical contact
- 1 remote control
CAUTION ! the remote control is programmed at delivery and ready for use. Don't take account of the instructions on the label stuck on it.


4 Display and controls

4.1 Device

- 1 Green On / Off button and for selection and save of the configuration
- 2 Battery housing
- 3 Remote control ; only the keys 0 to 9 are used
- 4 LCD display : displaying of elapsed time, number of shocks simulated, battery autonomy, text messages
- 5 Orange button for shock simulation according to the configuration of the device
- 6 Infrared receiver
- 7 Green diode for correct operation simulation
- 8 Orange diode for electrode fault simulation and scenario selection by remote control
- 9 Connector for connecting cable of the training electrodes

4.2 Electrodes and cable

- 10 Adhesive training electrodes; single use
CAUTION! : the electrodes are not designed to connect a patient electrically to the device; they are isolated. Don't use them with the defibrillator FRED easy!
- 11 Connector of the cable with a magnet to recognise the electrodes
- 12 Connection cable between electrodes and device
CAUTION! : The cable doesn't be used with the defibrillator FRED easy!
- 13 Discs of the cable to stick on the electrodes
- 14 Self adhesive discs of the electrodes

5 Scenario selection procedure

5.1 List of scenarios

SNR	Shock not recommended
SR	Shock recommended
Elect prob	Electrode problem

Scenario 0	SNR
Scenario 1	SNR / SR
Scenario 2	SNR / SR / SNR
Scenario 3	SR
Scenario 4	SR / SNR
Scenario 5	SR / SNR / SR
Scenario 6	Elect prob / SR
Scenario 7	SR/ Elect prob / SR
Scenario 8	SR / SR / SNR
Scenario 9	Infrared control function

Note : the last step of each scenario is indefinitely repeated.

5.2 Selection of a scenario

5.2.1 Selection of a scenario on the device

- Press the green button (1) to switch the device on. It functions with the last selected scenario.
- Place the device on standby mode by pressing the green button (1) for a few seconds. The green diode (7) flashes.
- Press at first the orange button (5), keep it pressed and press then the green button (1) to access in the configuration menu.

The selection in a menu is made by the green button; the passage to the next stage is made by the orange button.

- Line 1 : choice between automatic and semiautomatic functioning of the device
- Line 2 (only in the semiautomatic mode) : choice between automatic and manual analysis
- Line 3 : selection of one of the 9 scenarios (0 to 8) or selection of the infrared remote control.
- Line 4 : saving of the selected configuration and exit (with the green button)

- Select in this way the scenario according to the configuration of the FRED easy which you wish to use for training. Validate your selection by pressing the green button (1). The device then switches off and goes in standby mode. The green diode (7) flashes to inform the user that the device is operating normally.

The device can then be switched on normally and it will operate according to the selected scenario.

5.2.2 Selection of a scenario with the infrared remote control

The trainer has the possibility to select a scenario with the remote control when the device is in the standby mode (device off and green diode is flashing).

The orange diode (8) under the electrode connector flashes in the moment of the scenario selection to inform the user that the scenario has been taken in account by the device.

The device can then be switched on normally and it will operate according to the selected scenario.

In the case of the selection of scenario number 9, the device waits for the remote control's orders to change its status. If the device doesn't receive orders from the remote control the shock will always be recommended like in scenario 3.

5.3 Description of the functions associated with keys on the infrared remote control

Key 0	Immediate move to the next step
Key 1	Restart at the analysis step; reset of the counters of time and number of delivered shocks
Key 2	Shock recommended
Key 3	Shock not recommended
Key 4	The rhythm has changed
Key 5	Press on the electrodes (bad contact)
Key 6	Simulation of device off and green diode (7) off CAUTION ! The device is not really off. Press on one of the keys 0 to 5 to continue the training.
Key 7	Simulation of device off and green diode (7) is flashing CAUTION ! The device is not really in the standby mode. Press on one of the keys 0 to 5 to continue the training.
Key 8	Resume CPR
Key 9	Pause : the current scenario is interrupted : the display is frozen and the green diode (7) does not flash. Press again on key 9 to continue the training

5.4 To switch the device off

Press at the same time for a few seconds the orange button (5) and the green button (1), then slacken them to switch the device completely off. It is then no longer consuming energy.

CAUTION :
The device is only completely switched off when the green diode is no longer flashing. Do not leave it in standby mode, as this empties the battery unnecessarily.

6 Safety instructions

The following points shall be specially respected :

- The device has been designed for training purposes **only**. It may in any case not be used for defibrillation.
- The device is used with training electrodes.
- These instructions for use describe the complete FRED easy TRAINER device with all its optional features.
- The device is not designed for any use other than that described in these instructions. Such use is forbidden.
- The user shall be entirely responsible for use of other accessories than from SCHILLER.
- The user must be familiar with the functioning of the device.
- Before each use, the user must ensure the device is working correctly.
- The device is not designed for use in areas where an explosion hazard may occurs

7 Cleaning and maintenance

The outer surfaces of the device may be cleaned with cotton wool moistened with alcohol.

CAUTION !

**Risk of damage to the device
Never use phenol based cleaning agents or agents containing
peroxide derivatives to disinfect the surfaces of the device housing**

During cleaning, the batteries must be removed from the device. No liquid may be allowed to enter the device. If liquid does penetrate into the device, clean the device and have it verified completely.

8 Servicing and disposal at the end of the service life

Check the device performance once a year. At the end of its service life, the device and the accessories must be disposed of in compliance with the local regulations. Apart from the plug-in battery, the device does not contain hazardous material and can be disposed of like any other electronic equipment.

9 Specifications

- Operating mode : training in semi-automatic and automatic defibrillation
- Power supply :

Alkaline battery 6 Volts composed of 4 alkaline batteries LR6 1,5 V

CAUTION !
Do not recharge ! Do not put in the charger !
Do not use the battery in the defibrillator FRED easy !

Rechargeable Ni /Cd battery 6 Volts, 650 mAh

CAUTION !
Do not use the battery in the defibrillator FRED easy !


Battery charger 100-240 VAC, 50-60 Hz, 0,35A max ; Class II ; Charge time of an empty battery : 1h

Battery adapter

- 1 to 100% on-screen battery indicator corresponding to the battery voltage range from 4.5 Volts minimum to 6 Volts.

- Power consumption of the device : 35 mA without voice; approximately 250 mA with voice.
- Detection of electrode presence by means of a magnet placed in the connector of the connecting cable for the training electrodes.

9.1 Physical and environment requirements

- Transport /storage:

Temperature : – 30...+ 50°C

Relative humidity rate : 30...95 %, without condensation

Atmospheric pressure : 500...1060 hPa

- Operation:

Temperature : 0...+ 50 °C

Relative humidity rate : 30...95 %, without condensation

Atmospheric pressure : 700...1060 hPa

Dimensions and weight

Width 230 mm

Depth 220 mm

Height 70 mm

Weight < 1 kg

10 Warranty

SCHILLER guarantees its devices to the purchaser against all material or manufacturing defects for a period of one year as from the original purchase date. During this one-year period, SCHILLER, as soon as it receives a device presenting a material or manufacturing defect returned accompanied by written notification of the defect, will repair or replace the defective parts or the complete product if necessary.

All postal, shipping or handling costs will be borne exclusively by the user.

The manufacturer shall supply, on request, circuit diagrams, lists of components, descriptions, calibration instructions or any other information, for use by the qualified technical personnel of the user for repairing such parts of the device as have been stated to be repairable by the manufacturer. Such supply shall not in any case constitute permission or approval to modify or repair the devices.


All rights are reserved for the circuits, processes and names mentioned in these instructions for use.

SCHILLER shall not be held liable as a result of this warranty for any secondary or indirect damage, or in the event of unauthorised repairs or modifications having been made or undertaken, or when the whole or a part of the product has been damaged by accident, incorrect use or deliberate misuse. This warranty does not cover the battery, normal wear and tear, marks or any other irregularity of appearance which does not prevent or hinder the performance of the device.

11 Accessories

Part Number	Designation
6-02-0035	Case for 4 alkaline batteries (without batteries)
53108	Set of 4 alkaline batteries 1,5 Volt LR6
4-07-0010	Rechargeable NiCd 6 Volts battery
3-55-0030	Charger for 6 Volts NiCd battery, 100-240 VAC, 50-60 Hz Europe
0-21-0006	Reusable cable for training electrodes 0-21-0005. Length : 2 m.
0-21-0005	Adhesive training electrodes (10 sets of 2 electrodes), single use.
0-99-0001	Remote control, ready for use.

12 General algorithm

Note :

⦿ 1 ⦿ 20 seconds = Text 1 repeated every 20 seconds

Rem 1 : The step "PRESS ANALYSE" (text 2) is cancelled if the analyse is automatic or if the device functions in the automatic mode.

Rem 2 : The text 14 replaces the text 5 if the device functions in the automatic mode.

13 Spoken and displayed texts

	Spoken text	Text displayed
1	Place electrodes on chest and plug into machine.	PLACE ON ELECTRODES PLUG INTO MACHINE
2	Press green button.	PRESS GREEN BUTTON
3	Do not touch the patient, analysing.	DO NOT TOUCH THE PATIENT ANALYSING
4	Do not touch the patient, charging.	DO NOT TOUCH THE PATIENT CHARGING
5	Stand clear of patient, press orange button.	PRESS ORANGE BUTTON TO SHOCK
6	No shock advised	NO SHOCK ADVISED
7	Heart rhythm has changed.	RHYTHM HAS CHANGED
8	Immediately resume CPR. 30 chest compressions then 2 rescue breaths; continue until patient is breathing normally.	30 CHEST COMPRESSIONS THEN 2 RESCUE BREATHS
9	Movement detected. Stand clear.	MOVEMENT DETECTED STAND CLEAR
10	Check the electrodes.	CHECK THE ELECTRODES
14	Warning! Shock will be delivered 5-4-3-2-1.	WARNING ! SHOCK

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